

3M™ Littmann® Special Edition Stethoscopes

NEW! Rainbow Finish



Distinctive Finishes for a Distinctive You

Copper/Chocolate



Brass/Black



Gold/Black



Black/Black



Precision and color
combine in a unique
line of special edition
stethoscopes

**For those who have earned the right
to wear the best.**

You didn't choose the easiest career – you chose one of the most important. And by working hard every step of the way, you've made a difference.

It's time to mark your achievement – with a 3M™ Littmann® Special Edition Stethoscope.

Now choose from new rainbow or smoke finish chestpieces – or from any of the colors in the special edition line – and make your new Littmann stethoscope distinctively your own.


New smoke finish chestpiece. See back page.

**Precision.
Performance.
And a touch
of flair!**

New 3M™ Littmann® Special Edition Stethoscopes combine craftsmanship and style in a unique line of world-class instruments engineered to the exacting standards of medical professionals like you.

From the trusted leader in innovative auscultation technology, the Littmann brand name is your assurance of precision, acoustical superiority, innovative design and exceptional performance. Whichever model you choose, you can depend on Littmann stethoscopes to provide years of superior performance and service.

Make the choice that's right for you.



The quality of the instruments you use everyday sends a subtle but unmistakable message about your high professional standards.

To order, see your
3M Littmann
authorized dealer.

www.littmann.com



**3M Health Care
Littmann Brand Products**
3M Center, Building 275-4E-01
St. Paul, MN 55144-1000
USA
1-800 228 3957
www.3M.com/healthcare
www.littmann.com

Cat #	Product Description	Size Inches (cm)	
3M™ Littmann® Master Cardiology™ Stethoscopes			
2161	Black Finish Chestpiece, Black Tube	27 (69)	<ul style="list-style-type: none"> • Handcrafted chestpiece with the portable convenience of a single-side stethoscope • The tunable diaphragm conveniently alternates between low- and high-frequency sounds with a simple pressure change on the chestpiece • Dual-lumen design eliminates tube noise artifact from two tubes rubbing together • Seven-year warranty
NEW! 2175	Brass Finish Chestpiece, Black Tube	27 (69)	
NEW! 2176	Smoke Finish Chestpiece, Black Tube	27 (69)	
3M™ Littmann® Cardiology III™ Stethoscopes			
3128BRS	Brass Finish Chestpiece, Black Tube	27 (69)	<ul style="list-style-type: none"> • Two tunable diaphragms allow the user to conveniently alternate between the low- and high-frequency sounds without turning over the chestpiece • Dual-lumen design eliminates tube noise artifact from two tubes rubbing together • Five-year warranty
3137CPR	Copper Finish Chestpiece, Chocolate Tube	27 (69)	
3131BE	Black Finish Chestpiece, Black Tube	27 (69)	
3M™ Littmann® Master Classic II™ Stethoscopes			
2141	Black Finish Chestpiece, Black Tube	28 (71,1)	<ul style="list-style-type: none"> • Tunable diaphragm conveniently alternates between the low- and high-frequency sounds without turning over the chestpiece • Portable convenience of single-sided chestpiece • High-performance single-lumen tubing for excellent acoustic performance in a range of attractive colors • Three-year warranty
2142G	Gold Finish Chestpiece, Black Tube	28 (71,1)	
3M™ Littmann® Classic II S.E. Stethoscopes			
2201BRS	Brass Finish Chestpiece, Black Tube	28 (71,1)	<ul style="list-style-type: none"> • Features two-sided chestpiece that offers high acoustic sensitivity for superior performance • Tunable diaphragm on one side augments the performance of the traditional bell on the opposite side • Single-lumen tubing • Three-year warranty
2820CPR	Copper Finish Chestpiece, Chocolate Tube	28 (71,1)	
2218BE	Black Finish Chestpiece, Black Tube	28 (71,1)	
NEW! 2823	Rainbow Finish Chestpiece, Caribbean Blue Tube	28 (71,1)	
3M™ Littmann® Classic II Stethoscopes (Pediatric & Infant)			
NEW! 2153	(Pediatric) Rainbow Finish Chestpiece, Caribbean Blue Tube	28 (71,1)	<ul style="list-style-type: none"> • Traditional bell combined with a tunable diaphragm (on pediatric scopes only) that allows the user to conveniently alternate between low- and high-frequency sounds without turning over the chestpiece • Single-lumen tubing • Three-year warranty
NEW! 2157	(Infant) Rainbow Finish Chestpiece, Raspberry Tube	28 (71,1)	

**3M Health Care
3M Canada**
Post Office Box 5757
London, Ontario N6A 4T1
Canada
1-800 364 3577

Please recycle. Printed in U.S.A.
© 3M 2010. All rights reserved.
Issued 10/10. 7539HB
70-2010-7471-6

3M, Littmann, the Littmann logo, Master Cardiology and the shape of the Master Cardiology stethoscope chestpiece, Cardiology III and Master Classic II are trademarks of 3M.